

Analytical Chemistry

Mass Spectrometry

Prepared by Chan Man Hoong, edited by Lillian Lee

Instrument





Figure 2.1 – diagram of a mass spectrometer



Mass Spectrometer

- Six processes occur in a mass spectrometer:
- 1. Vapourisation chamber: If not already a gas, the compound is **vaporised in an oven.**
- 2. Ionisation chamber: **Electrons are fired at the gaseous molecules.** These **knock off other electrons from some of the molecules.**



Mass Spectrometer

- 3. Acceleration chamber: The gaseous ions are **accelerated by passing through an electric field**.
- 4. Electrostatic analyser: select ions of kinetic energy within a narrow range by using an electric field.
- 5. Deflection chamber: The fast-moving ions now pass through the poles of an electromagnet, where **they are deflected**.
- 6. The deflected ions **pass through a narrow slit and are collected on a metallic plate connected to an amplifier**.

Determination Using Mass Spectrometry

- Is used to identify unknown or new compounds.
- When a molecule is ionised it forms a **MOLECULAR ION** which can also undergo **FRAGMENTATION** or **RE-ARRANGEMENT** to produce particles of smaller mass.
- Only particles with a positive charge will be deflected and detected.
- The resulting spectrum has many peaks.
- The peak at **highest m/e ratio** is caused by **molecular ion** and indicates the **molecular mass**.
- The rest of the spectrum provides information about the structure.

Mass spectrum

- A mass spectrum is produced, which plots **relative abundance** against **mass/charge (m/e) ratio**.
- **Most ions that are formed in a mass spectrometer have a charge of +1** → x-axis is a measure of the mass.
- The **base peak** → from *particularly stable fragment (most abundant)* of the molecule.

Uses of Mass Spectrometry

Determination of the structures of organic compounds:

1. Measuring the relative heights of the molecular ion (**M**) peak and the (**M+1**) peak → can determine the **number of carbon atoms** in a molecule.
2. Using the (**M+2**) and (**M+4**) peaks → can identify halogen-containing (**Cl or Br**) compounds.
3. By measuring the **accurate mass** of a molecular ion → can determine its **molecular formula**.
4. By identifying the **fragments** produced when an ion breaks up inside a mass spectrometer → can often piece together the **structure** of the parent molecule.

Use of (M + 1) Peak

- Use to determine number of carbon atoms in molecule.
- Naturally-occurring carbon is composed of 98.9% ^{12}C and 1.1% ^{13}C (along with extremely small, and variable, amount of ^{14}C).
- The formula relating the (M+1)/(M) ratio the number of carbon atoms is:

$$n = \frac{100}{1.1} \left(\frac{A_{M+1}}{A_M} \right)$$

where n = number of carbon atoms
 A_{M+1} = the abundance of the M+1 peak
 A_M = the abundance of the molecular ion, M, peak.

Example 1(pg 58)

Compound A has a molecular ion at an m/e value of 120, and relative abundance 23%, and a peak at m/e 121 with a relative abundance of 2%. How many carbon atoms are in a molecule of A?

SAQ 1. pg 58

- Compound B contains carbon, hydrogen and oxygen only. Its mass spectrum contains a molecular ion peak at $m/e = 102$ (relative abundance 35%) and an $M+1$ peak at $m/e 103$ (1.5%).

Calculate the number of carbon atoms in the molecule, and hence deduce the number of oxygen atoms it contains, and its molecular formula.


Fragmentation patterns

- In ionisation chamber, bombardment by high speed *electrons* produces positive ions.
- The molecular ions formed undergo bond (homolytic) fission to produce molecular fragments:
 - 1 fragment is positive ,
 - 1 is a radical.
- $M \rightarrow M^+ \rightarrow B\bullet + A^+$

A^+ appear as further peaks.

$B\bullet$ does not appear as peaks in the spectrum.

e.g. Propanone



Fragmentation patterns

- Can be used to identify the structural formulae of compounds especially those with the same molecular formula.
- E.g. Propanone CH_3COCH_3 and
Propanal $\text{CH}_3\text{CH}_2\text{CHO}$
- Molecular of $\text{C}_3\text{H}_6\text{O}$ with $m/e = 58$

Mass Spectrum of propanone (pg 60)




Figure 2.2 – mass spectrum of propanone

Mass Spectrum of propanal (pg 60)




Figure 2.3 – mass spectrum of propanal

Mass Spectrum of propanone

- Highest $m/e = 58$
due to $\text{CH}_3\text{COCH}_3^+$ (whole molecule)
- $m/e = 43$ (base peak)
due to $(M - 15 = M - \text{CH}_3)$
i.e. CH_3CO^+
- $m/e = 15$ due to CH_3^+

(M), (M+2) and (M+4) peaks


- (M+2) and (M+4) peak occurs when compound contains Cl or Br atom.
- Both Cl & Br naturally occur as mixtures of 2 isotopes, with the relative abundances shown below:

element	isotope	relative abundance	approximate ratio
chlorine	^{35}Cl	75.8%	3:1
	^{37}Cl	24.2%	
bromine	^{79}Br	50.5%	1:1
	^{81}Br	49.5%	


(M+2) Peak

The **ratio of the M:(M+2) peak** should **reflect the natural abundances** given in the table (i.e. 3:1 for chlorine; 1:1 for bromine)


Because the two isotopes are of similar abundance, the peaks are of similar height.



- Example: mass spectrum of 1-bromopropane
- Note the isotope pattern at 122 and 124 that represent the M and M+2 in a 1:1 ratio.
- Loss of ^{79}Br from 122 or ^{81}Br from 124 gives the base peak a $m/z = 43$, corresponding to the propyl cation.
- Note that other peaks, such as those at $m/z = 107$ and 109 still contain Br and therefore also show the 1:1 isotope pattern.



- The first MS is of 2-chloropropane.
- Note the isotope pattern at 78 and 80 that represent the M and M+2 in a 3:1 ratio.
- Loss of ^{35}Cl from 78 or ^{37}Cl from 80 gives the base peak a $m/z = 43$, corresponding to the secondary propyl cation.
- Note that the peaks at $m/z = 63$ and 65 still contain Cl and therefore also show the 3:1 isotope pattern.



M+4 peaks

- If the molecule contains **two** chlorine atoms, (or two bromine atoms, or one of each) we should expect to see **three** molecular ions, at m/e values of **M**, **M+2** and **M+4**.
- E.g. M ($^{79}\text{Br}_2^+$), $M+2$ ($^{79}\text{Br}^{81}\text{Br}^+$) and $M+4$ ($^{81}\text{Br}_2^+$).
- The abundance ratio of Br_2 is 1:2:1
- The abundance ratio of Cl_2 is 9:6:1

Example (pg 58)

- dibromomethane, CH_2Br_2 ,
- three molecular ion peaks at m/e 172, 174 and 176 in the ratio 1:2:1,

Example (pg 58)

- dichloromethane, CH_2Cl_2 ,
- ratio 9:6:1,
- due to the greater natural abundance of the ^{35}Cl isotope.

SAQ 2 (pg 58)

- Calculate the $M : M+2 : M+4$ ratio for CH_2BrCl

High resolution mass spectra

- Can distinguish between ions that have same mass on low-resolution mass spectrum.

E.g :

name	structure	molecular formula
pentene	$\text{CH}_3\text{CH}_2\text{CH}_2\text{CH}=\text{CH}_2$	C_5H_{10}
aminopropanonitrile	$\text{CH}_3\text{CH}(\text{NH}_2)\text{CN}$	$\text{C}_3\text{H}_6\text{N}_2$
but-1-ene-3-one	$\text{CH}_2=\text{CHCOCH}_3$	$\text{C}_4\text{H}_6\text{O}$

- All have an approximate M_r of 70:

Molecular formulae from accurate masses

element	accurate relative atomic mass
H	1.0078
C	12.000
N	14.003
O	15.995

- These accurate isotopic mass can accurately measure the mass of molecular ion → only correspond to one molecular formula.

Accurate Masses

- $\text{C}_5\text{H}_{10} = 5 \times 12.000 + 10 \times 1.0078$
 $= 70.078$

- $\text{C}_3\text{H}_6\text{N}_2 = 3 \times 12.000 + 6 \times 1.0078$
 $+ 2 \times 14.003 = 70.053$

- $\text{C}_4\text{H}_6\text{O} = 4 \times 12.000 + 6 \times 1.0078$
 $+ 15.995 = 70.042$

- The last two differ by about 0.13%. This is well within the capabilities of a high-resolution mass spectrometer.

SAQ3 (pg 59)


- Explain whether a molecule having an accurate mass of 60.0574 is 1,2-diaminoethane, $C_2H_8N_2$, or propan-1-ol, C_3H_8O .

SAQ 4 (pg 61)


- Use the values of accurate relative atomic masses in the table to see whether it would be possible to decide whether the peak at $m/e = 29$ is due to $CH_3CH_2^+$ or CHO^+ .
- $CH_3CH_2^+ = (2 \times 12.000) + (5 \times 1.0078) = 29.039$
- $CHO^+ = 12.000 + 1.0078 + 15.995 = 29.003$

Molecules with same molecular formula but different structure

simplified mass spectrum of pentan-2-one - $\text{CH}_3\text{COCH}_2\text{CH}_2\text{CH}_3$



simplified mass spectrum of pentan-3-one - $\text{CH}_3\text{CH}_2\text{COCH}_2\text{CH}_3$



The molecular formula of B is $\text{C}_8\text{H}_8\text{O}_2$





Figure 2.4 – mass spectrum of compound B

Compound B


- **M:M+1 ratio of 11.5 : 1 suggests that it contains 8 carbon atoms,**
- **Accurate determination of its relative atomic mass suggests its molecular formula is $\text{C}_8\text{H}_8\text{O}_2$.**
- B could be *phenyl ethanoate, methyl benzoate or methyl benzoic acid*

Compound B


- phenyl ethanoate



- methyl benzoate




- methyl benzoic





Ease of formation of carbonium ion fragments

- $\text{CH}_3^+ < \text{RCH}_2^+ < \text{R}_2\text{CH}^+ < \text{R}_3\text{C}^+$



IDENTIFY THE COMPOUND



Common Small Ions		Common Neutral Fragments	
m/z	composition	mass loss	composition
15 amu	CH ₃	1 amu	H
17	OH	15	CH ₃
18	H ₂ O	17	OH
19	H ₂ O, F	18	H ₂ O
26	C ₂ H ₂ , CN	19	F
27	C ₂ H ₃	20	HF
28	C ₂ H ₄ , CO, H ₂ CN	27	C ₂ H ₃ , HCN
29	C ₂ H ₅ , CHO	28	C ₂ H ₄ , CO
30	CH ₂ NH ₂	30	CH ₂ O
31	CH ₃ O	31	CH ₃ O
33	SH, CH ₂ F	32	CH ₃ O, S
34	H ₂ S	33	CH ₃ + H ₂ O, HS
35(37)	Cl	33	H ₂ S
36(38)	HCl	35(37)	Cl
39	C ₃ H ₃	36(38)	HCl
41	C ₃ H ₅ , C ₂ H ₅ N	42	C ₃ H ₅ , C ₂ H ₅ O, C ₂ H ₅ N
42	C ₃ H ₆ , C ₂ H ₆ O, C ₂ H ₅ N	43	C ₃ H ₇ , CH ₃ CO
43	C ₃ H ₇ , CH ₃ CO	44	CO ₂ O, CONH ₂
44	C ₂ H ₄ O	45	C ₂ H ₄ O
46	NO ₂	55	C ₄ H ₇
56	C ₄ H ₈	57	C ₄ H ₉
57	C ₄ H ₉	59	C ₄ H ₉ O ₂
60	CH ₃ CO ₂	60	C ₂ H ₅ O ₂
79(81)	Br	64	SO ₂
80(82)	HBr	79(81)	Br
91	C ₄ H ₇	80(82)	HBr
127	I	127	I
128	HI	128	HI

APPLICATIONS

- By coupling a GLC in conjunction with a mass spectrometer, **rapid analysis of complex mixtures** is possible.
- This is particularly useful for **determining the products and relative yields from organic reactions** and for **monitoring industrial processes**.

- ◆ Mass spec. are **suitable for analyzing volatile compounds**.
- ◆ It can also be used to **analyze proteins and polypeptides**. This is achieved by methylating the -N-H groups which disrupts the HB, hence more volatile.
- ◆ This allows a very rapid method of **determining the a. a sequence in the molecules**.
- ◆ This technique is usually computer-linked and has the advantage that very small quantities are required and sequences of amino acids may be rapidly established.